


Phospho-Ezh2/Kmt6 (T311) Rabbit Polyclonal Antibody For Wb, Elisa - Q15910

Item Number: CM0001206



Introduction

Phospho-specific rabbit polyclonal antibody detecting EZH2 phosphorylated at Thr311. Validated for WB and ELISA, suitable for epigenetic, chromatin, and cancer research. Cross-reacts with rat.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EZH2/KMT6-T311 Rabbit pAb
Remarks/Alias	WVS; ENX1; KMT6; WVS2; ENX-1; EZH2b; KMT6A; Phospho-EZH2/KMT6-T311
Species	Human
GeneID (human)	2146
GeneID	2146
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around T311 of human EZH2/KMT6 (NP_001190176.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q15910
Protein Weight	85kDa
Shipping	Ice bag