


Phospho-Ezr-T567/Rdx-T564/Msn-T558 Rabbit Pab For Wb, Elisa - P15311 / P35241 / P26038

Item Number: CM0001008



Introduction

Phospho-specific rabbit polyclonal antibody detecting Ezrin (T567), Radixin (T564), and Moesin (T558). Validated for WB and ELISA in human and rat. Useful for studying ERM protein activation and cytoskeleton-membrane interactions. Predicted MW 23-71 kDa.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EZR-T567/RDX-T564/MSN-T558 Rabbit pAb
Species	Human
GeneID (Human)	7430/4478/5962
GeneID	7430/4478/5962
Immunogen	A synthetic phosphorylated peptide around T567 of human EZR/RDX/MSNEZR (NP_003370.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P15311/P35241/P26038
Protein Weight	23-71 kDa
Shipping	Ice bag