


Phospho-Ezrin-T567 Rabbit Pab For Wb - P15311

Item Number: CM0001107



Introduction

Phospho-specific rabbit polyclonal antibody targeting human Ezrin at threonine 567. Validated for Western blot (WB) in human, mouse, and rat. Suitable for investigating cytoskeletal dynamics, cell adhesion, and membrane-cytoskeleton interactions.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Ezrin-T567 Rabbit pAb
Remarks/Alias	CVL; CVIL; VIL2; HEL-S-105; Phospho-Ezrin-T567
Species	Human
GeneID (Human)	7430
GeneID	7430
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding T567 of human Ezrin
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P15311
Protein Weight	69kDa
Shipping	Ice bag