


Phospho-Ezrin (T567) + Radixin (T564) + Moesin (T558) Rabbit Monoclonal Antibody For Wb And Elisa - P15311 / P35241 / P26038

Item Number: CM0001332



Introduction

Phospho-specific rabbit monoclonal antibody detecting phosphorylated Ezrin (T567), Radixin (T564), and Moesin (T558) in human, mouse, rat. Suitable for WB and ELISA. Ideal for studying ERM protein activation and cytoskeletal signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Ezrin-T567+Radixin-T564+Moesin-T558 Rabbit mAb
Remarks/Alias	CVL; CVIL; VIL2; HEL-S-105; DFNB24; Radixin; HEL70; IMD50; Moesin
Species	Human
GenelD (Human)	7430/5962/4478
GenelD	7430/5962/4478
Immunogen	A synthetic phosphorylated peptide around T567/T564/T558 of human Ezrin/Radixin/Moesin
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15311/P35241/P26038
Protein Weight	23-71 kDa
Shipping	Ice bag