


Phospho-Ezrin-Thr567/Radixin-Thr564/Moesin-Thr558 Rabbit Pab For Wb, Elisa - P15311 / P35241 / P26038

Item Number: CM0001006



Introduction

Phospho-specific rabbit polyclonal antibody detecting Ezrin (Thr567), Radixin (Thr564), and Moesin (Thr558) for WB and ELISA assays in human samples. Ideal for studying cytoskeleton-membrane interactions and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Ezrin-Thr567/Radixin-Thr564/Moesin-Thr558 Rabbit pAb
Remarks/Alias	CVL; CVIL; VIL2; HEL-S-105; DFNB24; HEL70; IMD50; Phospho-Ezrin-Thr567/Radixin-Thr564/Moesin-Thr558
Species	Human
GeneID (Human)	7430/5962/4478
Immunogen	A synthetic phosphorylated peptide around T567 of human Ezrin/Radixin/Moesin (NP_003370.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15311/P35241/P26038
Protein Weight	23-71 kDa
Shipping	Ice bag