


Phospho-Ezrin-Y353 Rabbit Pab For Wb And Elisa - P15311

Item Number: CM0001093



Introduction

Phospho-specific rabbit polyclonal antibody targeting human Ezrin phosphorylated at Y353 (SWISS P15311), validated for Western blot and ELISA applications, ideal for cytoskeleton and cell polarity studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Ezrin-Y353 Rabbit pAb
Remarks/Alias	CVL; CVIL; VIL2; HEL-S-105; Phospho-Ezrin-Y353
Species	Human
GeneID (human)	7430
GeneID	7430
Immunogen	A synthetic phosphorylated peptide around Y353 of human Ezrin (NP_003370.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15311
Protein Weight	69kDa
Shipping	Ice bag