


Phospho-Fadd-S194 Rabbit Pab For Wb, Elisa - Q13158

Item Number: CM0000516



Introduction

Phospho-specific rabbit polyclonal antibody targeting human FADD phosphorylated at S194, validated for Western blot and ELISA. Ideal for apoptosis and innate immunity research applications.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FADD-S194 Rabbit pAb
Remarks/Alias	GIG3; IMD90; MORT1; Phospho-FADD-S194
Species	Human
GeneID (human)	8772
GeneID	8772
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S194 of human FADD (NP_003815.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13158
Protein Weight	23kDa
Shipping	Ice bag