

Phospho-Fak-Y397 Rabbit Mab For Wb And Elisa - Q05397

Item Number: CM0000863



Introduction

Phospho-specific rabbit monoclonal antibody against FAK phosphorylated at Y397, validated for Western blot and ELISA applications in human samples. Ideal for studying integrin signaling, cell migration, and focal adhesion dynamics.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FAK-Y397 Rabbit mAb
Remarks/Alias	FAK; FADK; FAK1; FRNK; FADK 1; PPP1R71; p125FAK; pp125FAK; Phospho-FAK-Y397
Species	Human
Gene ID	5747
Immunogen	A synthetic phosphorylated peptide around Y397 of human FAK (NP_722560.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q05397
Protein Weight	119kDa
Shipping	Ice bag