


Phospho-Fak-Y397 Rabbit Polyclonal Antibody For Wb, Elisa - Q05397

Item Number: CM0001945



Introduction

Phospho-specific rabbit polyclonal antibody targeting FAK phosphorylated at Tyr397. Validated for WB and ELISA in human, mouse, rat samples. Ideal for studying integrin signaling and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FAK-Y397 Rabbit pAb
Remarks/Alias	FAK; FADK; FAK1; FRNK; FADK 1; PPP1R71; p125FAK; pp125FAK; Phospho-FAK-Y397
Species	Human
GeneID (human)	5747
GeneID	5747
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05397
Protein Weight	119kDa
Shipping	Ice bag