


Phospho-Fak-Y576/577 Rabbit Polyclonal Antibody For Wb, Elisa - Q05397

Item Number: CM0001959



Introduction

Phospho-specific rabbit polyclonal antibody targeting FAK (phosphorylated at Tyr576/577) for Western blot and ELISA applications in human samples. Useful for studying integrin signaling, cell migration, and adhesion.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FAK-Y576/577 Rabbit pAb
Remarks/Alias	FAK; FADK; FAK1; FRNK; FADK 1; PPP1R71; p125FAK; pp125FAK; Phospho-FAK-Y576/577
Species	Human
Gene ID (Human)	5747
Gene ID	5747
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around Y576 & Y577 of human PTK2 (NP_722560.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q05397
Protein Weight	119kDa
Shipping	Ice bag