


Phospho-Fgfr3-Y724 Rabbit Pab - P22607

Item Number: CM0001115



Introduction

Rabbit polyclonal antibody specific for phosphorylated FGFR3 (Tyr724). Validated for WB and ELISA in human and mouse samples. Ideal for studying FGFR3 signaling in skeletal development and cancer.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FGFR3-Y724 Rabbit pAb
Remarks/Alias	ACH; CEK2; JTK4; CD333; HSGFR3EX; Phospho-FGFR3-Y724
Species	Human
GeneID (Human)	2261
GeneID	2261
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P22607
Protein Weight	88kDa
Shipping	Ice bag