

Phospho-Filamin A-S2152 Rabbit Pab For Wb, Elisa - P21333

Item Number: CM0001772



Introduction

Phospho-Filamin A-S2152 rabbit polyclonal antibody validated for WB and ELISA against human, mouse, rat Filamin A (281kDa). Ideal for actin cytoskeleton, signaling, and disease research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Filamin A-S2152 Rabbit pAb
Remarks/Alias	FLN; FMD; MNS; OPD; ABPX; CSBS; CVD1; FGS2; FLN1; NHBP; OPD1; OPD2; XLVD; XMVD; FLN-A; ABP-280; Phospho-Filamin A-S2152
Species	Human
Gene ID (Human)	2316
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21333
Protein Weight	281kDa
Shipping	Ice bag