

Phospho-Foxo1-S319 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12778

Item Number: CM0000813



Introduction

Phospho-specific rabbit polyclonal antibody targeting phosphorylation at Ser319 of human FOXO1 (Gene ID: 2308). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying insulin signaling and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FOXO1-S319 Rabbit pAb
Remarks/Alias	FKH1; FKHR; FOXO1A; Phospho-FOXO1-S319
Species	Human
GeneID (human)	2308
GeneID	2308
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S319 of human FOXO1 (NP_002006.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12778
Protein Weight	70kDa
Shipping	Ice bag