

Phospho-Foxo3A-S253 Rabbit Pab For Wb, Ihc-P, Elisa - O43524

Item Number: CM0000803



Introduction

Phospho-specific rabbit polyclonal antibody against FOXO3A (phospho-Ser253). Validated for WB, IHC-P, ELISA. Detects human, mouse, rat. Ideal for apoptosis, autophagy, and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FOXO3A-S253 Rabbit pAb
Remarks/Alias	FOXO2; AF6q21; FKHL1; FOXO3A; FKHL1P2; Phospho-FOXO3A-S253
Species	Human
GeneID (Human)	2309
GeneID	2309
Immunogen	Synthetic phosphorylated peptide around S253 of human FOXO3A
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43524
Protein Weight	71kDa
Shipping	Ice bag