

## Phospho-Foxo3A-S425 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43524

Item Number: CM0001987



### Introduction

Phospho-specific rabbit polyclonal antibody targeting FOXO3A phosphorylation at Ser425. Validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Supports apoptosis and autophagy research.

[Learn More](#)

| Parameter        | Value                                                                                           |
|------------------|-------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-FOXO3A-S425 Rabbit pAb                                                                  |
| Remarks/Alias    | FOXO2; AF6q21; FKHL1; FOXO3A; FKHL1P2; Phospho-FOXO3A-S425                                      |
| Species          | Human                                                                                           |
| GeneID (Human)   | 2309                                                                                            |
| GeneID           | 2309                                                                                            |
| Immunogen        | Synthetic peptide A synthetic phosphorylated peptide around S425 of human FOXO3A (NP_001446.1). |
| Source           | Rabbit                                                                                          |
| Category         | Phosphorylated Antibodies                                                                       |
| Application      | WB, IHC-P, ELISA                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                               |
| SWISS            | O43524                                                                                          |
| Protein Weight   | 71kDa                                                                                           |
| Shipping         | Ice bag                                                                                         |