

Phospho-Foxo4-S197 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P98177

Item Number: CM0000925



Introduction

Phospho-FOXO4 (S197) Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat FOXO4 phosphorylated at Ser197. Ideal for insulin signaling, cell cycle, and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FOXO4-S197 Rabbit pAb
Remarks/Alias	AFX; AFX1; MLLT7; Phospho-FOXO4-S197
Species	Human
GeneID (Human)	4303
GeneID	4303
Immunogen	A synthetic phosphorylated peptide around S197 of human FOXO4 (NP_005929.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98177
Protein Weight	54kDa
Shipping	Ice bag