


Phospho-Gata1-S310 Rabbit Polyclonal Antibody For Wb, Elisa - P15976

Item Number: CM0001605



Introduction

Phospho-GATA1-S310 rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat GATA1 phosphorylated at Ser310. Suitable for erythroid development and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GATA1-S310 Rabbit pAb
Remarks/Alias	GF1, GF-1, NFE1, XLTT, ERYF1, NF-E1, XLANP, XLTDA, GATA-1, HAEADA, Phospho-GATA1-S310
Species	Human
GeneID (Human)	2623
GeneID	2623
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding S310 of human GATA1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15976
Protein Weight	43kDa
Shipping	Ice bag