


Phospho-GluR1/Gria1 (S836) Rabbit Polyclonal Antibody For Wb - P42261

Item Number: CM0001100



Introduction

Phospho-specific rabbit polyclonal antibody targeting GluR1/GRIA1 at serine 836. Designed for Western blot detection of human, mouse, and rat AMPA receptor subunit. Ideal for neuroscience and synaptic plasticity research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GluR1/GRIA1-S836 Rabbit pAb
Remarks/Alias	GLUH1; GLUR1; GLURA; GluA1; HBGR1; MRD67; MRT76; Phospho-GluR1/GRIA1-S836
Species	Human
GeneID (Human)	2890
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding S836 of human GluR1/GRIA1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P42261
Protein Weight	102kDa
Shipping	Ice bag