


Phospho-GluR2/Gria2-S880 Rabbit Pab For Wb - P42262

Item Number: CM0001796



Introduction

Phospho-specific rabbit polyclonal antibody targeting GluR2/GRIA2 at Ser880. Validated for WB with human, mouse, rat cross-reactivity. Suitable for neuroscience research on AMPA receptor signaling and synaptic plasticity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GluR2/GRIA2-S880 Rabbit pAb
Remarks/Alias	GLUR2; GLURB; GluA2; HBGR2; NEDLIB; gluR-2; gluR-B; GluR-K2; Phospho-GluR2/GRIA2-S880
Species	Human
GeneID (Human)	2891
GeneID	2891
Immunogen	A phospho specific peptide corresponding to residues surrounding S880 of human GluR2/GRIA2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P42262
Protein Weight	99kDa
Shipping	Ice bag