

Phospho-Grin2B-Y1474 Rabbit Pab For Wb, Elisa - Q13224

Item Number: CM0000802



Introduction

Phospho-GRIN2B-Y1474 rabbit polyclonal antibody detects phosphorylated NMDA receptor subunit NR2B at Y1474 in human, mouse, and rat samples. Validated for WB and ELISA, this antibody supports neuroscience research of synaptic plasticity and excitotoxicity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GRIN2B-Y1474 Rabbit pAb
Remarks/Alias	NR3; MRD6; NR2B; hNR3; DEE27; EIEE27; GluN2B; NMDAR2B; Phospho-GRIN2B-Y1474
Species	Human
Gene ID (Human)	2904
Gene ID	2904
Immunogen	A synthetic phosphorylated peptide around Y1474 of human GRIN2B (NP_000825.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13224
Protein Weight	166kDa
Shipping	Ice bag