


Phospho-Gsk3B-S9 Rabbit Mab For Wb, Ihc-P, Elisa - P49841

Item Number: CM0001592



Introduction

Phospho-GSK3 β -S9 Rabbit mAb (P49841) for Western blot, IHC-P, ELISA. Detects human GSK3 β phosphorylated at Ser9. Ideal for studying Wnt signaling, insulin resistance, and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GSK3 β -S9 Rabbit mAb
Remarks/Alias	GSK3B; gsk-3 β ; Phospho-GSK3 β -S9
Species	Human
GeneID (human)	2932
GeneID	2932
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S9 of human GSK3 beta (P49841)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P49841
Protein Weight	47 kDa
Shipping	Ice bag