


Phospho-Gsk3B-S9 Rabbit Pab - P49841

Item Number: CM0000590



Introduction

Rabbit polyclonal antibody targeting phosphorylated Ser9 of human GSK3 β , validated for WB, IF/ICC, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for studying Wnt signaling and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GSK3 β -S9 Rabbit pAb
Remarks/Alias	GSK3B; gsk-3 β ; Phospho-GSK3 β -S9
Species	Human
Gene ID (Human)	2932
Gene ID	2932
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49841
Protein Weight	47kDa
Shipping	Ice bag