


Phospho-Gsk3B-Y216 + Gsk3A-Y279 Rabbit Monoclonal Antibody For Wb, Elisa - P49840/P49841

Item Number: CM0001346



Introduction

Phospho-specific rabbit monoclonal antibody detecting GSK3 β phosphorylated at Y216 and GSK3 α phosphorylated at Y279. Validated for Western blot and ELISA, showing cross-reactivity with mouse. Ideal for research in glucose metabolism, Wnt signaling, and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GSK3 β -Y216 + GSK3 α -Y279 Rabbit mAb
Remarks/Alias	GSK3B; GSK 3 beta; GSK3A; GSK 3 alpha; Phospho-GSK3 α -S21+GSK3 β -S9
Species	Human
GeneID (Human)	2931/2932
GeneID	2931/2932
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y216&Y279 of human GSK3 β /GSK3 α (NP_063937.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P49840/P49841
Protein Weight	51kDa/46kDa
Shipping	Ice bag