


Phospho-Gsk3B-Y216 + Gsk3A-Y279 Rabbit Pab - P49840 / P49841

Item Number: CM0001003



Introduction

Phospho-specific rabbit polyclonal antibody recognizing GSK3 β -Y216 and GSK3 α -Y279, suitable for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for studies on insulin signaling, Wnt pathway, and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GSK3 β -Y216 + GSK3 α -Y279 Rabbit pAb
Remarks/Alias	GSK3B; gsk-3 β ; Phospho-GSK3 β -Y216 + GSK3 α -Y279
Species	Human
Gene ID	2931 / 2932
Immunogen	Synthetic phosphorylated peptide around Y216/Y279 of human GSK3 β /GSK3 α
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49840 / P49841
Protein Weight	47kDa
Shipping	Ice bag