


## Phospho-Gsk3B-Y216 + Gsk3A-Y279 Rabbit Polyclonal Antibody For Wb, Elisa - P49841

Item Number: CM0001216



### Introduction

Phospho-specific antibody against GSK3 $\beta$  (pTyr216) and GSK3 $\alpha$  (pTyr279). Validated for WB and ELISA in human samples. Essential for Wnt signaling, glucose metabolism, and Alzheimer's disease research.

[Learn More](#)

| Parameter        | Value                                                                               |
|------------------|-------------------------------------------------------------------------------------|
| Product Name     | Phospho-GSK3 $\beta$ -Y216 + GSK3 $\alpha$ -Y279 Rabbit pAb                         |
| Remarks/Alias    | GSK3B; gsk-3 $\beta$ ; Phospho-GSK3 $\beta$ -Y216 + GSK3 $\alpha$ -Y279             |
| Species          | Human                                                                               |
| Gene ID (Human)  | 2931/2932                                                                           |
| Gene ID          | 2931/2932                                                                           |
| Immunogen        | A synthetic phosphorylated peptide around Y216 of human GSK3 $\beta$ (NP_001139628) |
| Source           | Rabbit                                                                              |
| Category         | Phosphorylated Antibodies                                                           |
| Application      | WB, ELISA                                                                           |
| Cross Reactivity | Human                                                                               |
| SWISS            | P49841                                                                              |
| Protein Weight   | 47kDa                                                                               |
| Shipping         | Ice bag                                                                             |