


Phospho-Hdac4-S246/Hdac5-S259/Hdac9-S220 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P56524 / Q9Uql6 / Q9Ukv0

Item Number: CM0000428



Introduction

Phospho-specific rabbit polyclonal antibody detecting HDAC4 (S246), HDAC5 (S259), and HDAC9 (S220) for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for chromatin regulation and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220 Rabbit pAb
Remarks/Alias	HDAC4/HDAC5/HDAC9; Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220
Species	Human
GenelD (human)	9759/10014/9734
GenelD	9759/10014/9734
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding S246 of human HDAC4
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56524/Q9UQL6/Q9UKV0
Protein Weight	57-65kDa/97-122kDa
Shipping	Ice bag