


Phospho-Hdac4-S632 Rabbit Monoclonal Antibody For Wb And Elisa - P56524

Item Number: CM0001615



Introduction

Phospho-HDAC4-S632 Rabbit mAb is a phosphorylation-specific antibody designed for Western blot (WB) and ELISA applications. It detects human HDAC4 phosphorylated at serine 632, supporting epigenetics and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HDAC4-S632 Rabbit mAb
Remarks/Alias	HD4; AHO3; BDMR; HDACA; HA6116; HDAC-4; HDAC-A; NEDCHF; NEDCHID; Phospho-HDAC4-S632
Species	Human
GeneID	9759
Immunogen	A synthetic phosphorylated peptide around S632 of human Phospho-HDAC4-S632 (NP_006028.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P56524
Protein Weight	119kDa
Shipping	Ice bag