


Phospho-Hdac4-S632 Rabbit Polyclonal Antibody For Wb - P56524

Item Number: CM0000786



Introduction

Phospho-HDAC4-S632 Rabbit Polyclonal Antibody detects endogenous levels of HDAC4 phosphorylated at serine 632. Validated for Western blot in human, mouse, and rat samples. Ideal for chromatin and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HDAC4-S632 Rabbit pAb
Remarks/Alias	HD4; AHO3; BDMR; HDACA; HA6116; HDAC-4; HDAC-A; NEDCHF; NEDCHID; Phospho-HDAC4-S632
Species	Human
GeneID (human)	9759
GeneID	9759
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P56524
Protein Weight	119kDa
Shipping	Ice bag