

Phospho-Hdac5-S498 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uql6

Item Number: CM0001466



Introduction

Rabbit polyclonal antibody against phosphorylated HDAC5 (S498) for Western blot and ELISA applications. Reacts with human, mouse, rat HDAC5. Suitable for research on epigenetic regulation, chromatin remodeling, and transcriptional repression.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HDAC5-S498 Rabbit pAb
Remarks/Alias	HD5, NY-CO-9, Phospho-HDAC5-S498
Species	Human
Gene ID (Human)	10014
Gene ID	10014
Immunogen	A synthetic phosphorylated peptide around S498 of human HDAC5 (NP_005465.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQL6
Protein Weight	122kDa
Shipping	Ice bag