


Phospho-Hdac8-S39 Rabbit Pab - Q9By41

Item Number: CM0000900



Introduction

Phospho-specific rabbit polyclonal antibody targeting HDAC8 at Ser39. Validated for IF/ICC on human samples; detects histone deacetylase 8 involved in chromatin remodeling and transcriptional repression.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HDAC8-S39 Rabbit pAb
Remarks/Alias	HD8; WTS; RPD3; CDA07; CDLS5; KDAC8; MRXS6; HDACL1; Phospho-HDAC8-S39
Species	Human
Gene ID (Human)	55869
Gene ID	55869
Immunogen	Synthetic peptide; A phospho specific peptide corresponding to residues surrounding S39 of human HDAC8
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IF/ICC
Cross Reactivity	Human
SWISS	Q9BY41
Protein Weight	42kDa
Shipping	Ice bag