


Phospho-Histone H2Ax-S139 Rabbit Mab For Wb, Ip, If/Icc, Elisa - P16104

Item Number: CM0000817



Introduction

Rabbit monoclonal antibody targeting phosphorylated Ser139 of Histone H2AX (H2A.X). Validated for WB, IP, IF/ICC, and ELISA with cross-reactivity in human, mouse, and rat. Ideal for DNA damage and repair research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H2AX-S139 Rabbit mAb
Remarks/Alias	H2A.X; H2A/X; H2AFX
Species	Human
Gene ID (Human)	3014
Gene ID	3014
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S139 of human Histone H2AX (NP_002096.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16104
Protein Weight	15 kDa
Shipping	Ice bag