


## Phospho-Histone H3-S10 Rabbit Mab For Wb, Ihc-P, Ip, Db, Elisa - Q16695 / P68431

Item Number: CM0000771



### Introduction

High-specificity Phospho-Histone H3 (Ser10) Rabbit Monoclonal Antibody validated for WB, IHC-P, IP, DB, ELISA. Reacts with human, mouse, rat histone H3; ideal for epigenetic and cell signaling studies.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-Histone H3-S10 Rabbit mAb                                                                                    |
| Remarks/Alias    | H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3 |
| Species          | Human                                                                                                                |
| GeneID (human)   | 8290/8350                                                                                                            |
| GeneID           | 8290/8350                                                                                                            |
| Immunogen        | Synthetic peptide                                                                                                    |
| Source           | Rabbit                                                                                                               |
| Category         | Phosphorylated Antibodies                                                                                            |
| Application      | WB, IHC-P, IP, DB, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                                                                      |
| SWISS            | Q16695/P68431                                                                                                        |
| Protein Weight   | 15 kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                              |