


Phospho-Histone H3-S10/T11 Rabbit Pab - Q16695 / P68431

Item Number: CM0001434



Introduction

Phospho-specific rabbit polyclonal antibody detecting Histone H3 phosphorylated at Ser10 and Thr11. Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat, and predicted wide range. Suitable for chromatin dynamics and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-S10/T11 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Phospho-Histone H3-S10/T11
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S10 & T11 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag