


Phospho-Histone H3 (Ser28) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16695 / P68431

Item Number: CM0001437



Introduction

Phospho-specific rabbit polyclonal antibody targeting Histone H3 phosphorylated at Ser28, suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat H3 with wide predicted cross-reactivity. Ideal for chromatin, cell cycle, and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-S28 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Phospho-Histone H3-S28
Species	Human
Gene ID (Human)	H3
Gene ID	8290/8350
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S28 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag