


Phospho-Histone H3-T11 Rabbit Pab - Q16695 / P68431

Item Number: CM0001030



Introduction

Phospho-Histone H3-T11 Rabbit pAb (Q16695/P68431) for WB, IHC-P, IF/ICC, ELISA. High specificity for human, mouse, rat, and wide-range predicted. Validated for detecting phosphorylated histone H3 at Thr11 in testicular and other cells.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-T11 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Phospho-Histone H3-T11
Species	Human
Gene ID (Human)	H3
Gene ID	8290/8350
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T11 of human HIST1H3B (NP_003520.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag