


Phospho-Histone H3-T3 Rabbit Pab For Wb, Elisa - Q16695 / P68431

Item Number: CM0000899



Introduction

Phospho-specific rabbit polyclonal antibody targeting Thr-3 phosphorylated Histone H3. Validated for WB and ELISA, it detects human, mouse, and rat H3 with predicted wide cross-reactivity. Ideal for chromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-T3 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Phospho-Histone H3-T3
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around T3 of human HIST1H3A (NP_003520.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag