


Phospho-Histone H4-S1 Rabbit Pab For Wb, Elisa - P62805

Item Number: CM0001131



Introduction

Phospho-Histone H4-S1 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human, mouse, rat, and other species. Suitable for chromatin and histone modification studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H4-S1 Rabbit pAb
Remarks/Alias	H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FA; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4A; Phospho-Histone H4-S1
Species	Human
GeneID (Human)	H4
GeneID	8359
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S1 of human Histone H4 (NP_003529.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag