


Phospho-Hsp27/Hspb1-S15 Rabbit Mab - P04792

Item Number: CM0001014



Introduction

Phospho-specific rabbit monoclonal antibody targeting HSP27/HSPB1 at Ser15. Validated for WB, IF/ICC, ELISA. Suitable for human samples, detects 23 kDa protein. Ideal for stress response and chaperone studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSP27/HSPB1-S15 Rabbit mAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; Phospho-HSP27/HSPB1-S15
Species	Human
GeneID (Human)	3315
GeneID	3315
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S15 of human HSP27/HSPB1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P04792
Protein Weight	23kDa
Shipping	Ice bag