


Phospho-Hsp27/Hspb1-S82 Rabbit Mab For Wb, Elisa - P04792

Item Number: CM0001766



Introduction

Phospho-specific rabbit monoclonal antibody against human HSP27/HSPB1 (Ser82). Validated for WB and ELISA; cross-reacts with human and rat. Ideal for stress response, chaperone function, and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSP27/HSPB1-S82 Rabbit mAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; Phospho-HSP27/HSPB1-S82
Species	Human
Gene ID	3315
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S82 of human HSP27/HSPB1 (P04792).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04792
Protein Weight	23 kDa
Shipping	Ice bag