


Phospho-Hsp27/Hspb1-S82 Rabbit Pab - P04792

Item Number: CM0000434



Introduction

Phospho-HSP27/HSPB1-S82 Rabbit pAb detects endogenous phospho-HSP27 at Ser82. Validated for Western blot, IHC-P, IP, and ELISA. Ideal for stress response, chaperone activity, and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSP27/HSPB1-S82 Rabbit pAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; Phospho-HSP27/HSPB1-S82
Species	Human
GeneID (Human)	3315
GeneID	3315
Immunogen	A synthetic phosphorylated peptide around S82 of human HSP27/HSPB1 (NP_001531.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	P04792
Protein Weight	23kDa
Shipping	Ice bag