

Phospho-Igf1R-Y1161 Rabbit Polyclonal Antibody For Wb And Elisa - P08069

Item Number: CM0000943



Introduction

Phospho-IGF1R-Y1161 Rabbit Polyclonal Antibody for WB and ELISA. Reactivity: human, mouse, rat. Targets phosphorylated Y1161 of IGF1R, a tyrosine kinase receptor regulating cell growth and tumor transformation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IGF1R-Y1161 Rabbit pAb
Remarks/Alias	IGFR; CD221; IGFR; JTK13; Phospho-IGF1R-Y1161
Species	Human
GeneID (Human)	3480
GeneID	3480
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y1161 of human IGF1R (NP_000866.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08069
Protein Weight	155kDa
Shipping	Ice bag