

Phospho-Igf1R-Y1165/1166/Insr-Y1189/1190 Rabbit Polyclonal Antibody - P08069 / P06213

Item Number: CM0000541



Introduction

Phospho-specific rabbit polyclonal antibody recognizing the activated forms of IGF1R (Y1165/Y1166) and INSR (Y1189/Y1190). Suitable for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for growth signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IGF1R-Y1165/1166/INSR-Y1189/1190 Rabbit pAb
Remarks/Alias	IGF1R/INSR; Phospho-IGF1R-Y1165/1166/INSR-Y1189/1190
Species	Human
Gene ID	3480 (IGF1R), 3643 (INSR)
Immunogen	A synthetic phosphorylated peptide around Y1165 & Y1166 of human IGF1R (NP_000866.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08069, P06213
Protein Weight	154kDa/155kDa/156kDa
Shipping	Ice bag