

Phospho-IkB α -S36 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P25963

Item Number: CM0001520



Introduction

Phospho-specific rabbit monoclonal antibody against I κ B α (Ser36). Suitable for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Validated for NF- κ B signaling pathway research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IkB α -S36 Rabbit mAb
Remarks/Alias	IKBA; MAD-3; NFKB1; EDAID2; Phospho-IkB α -S36
Species	Human
Gene ID (Human)	4792
Gene ID	4792
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S36 of human I κ B α (P25963).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25963
Protein Weight	36kDa
Shipping	Ice bag