

Phospho-Ikba-Y42 Rabbit Pab For Wb - P25963

Item Number: CM0001033



Introduction

Phospho-IkBa-Y42 Rabbit pAb specifically detects human NF-kappa-B inhibitor alpha phosphorylated at tyrosine 42 by Western blot. Ideal for studying NF-kB signaling, immune response, and inflammation.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-IkBa-Y42 Rabbit pAb
Remarks/Alias	IKBA; MAD-3; NFKBI; EDAID2; Phospho-IkBa-Y42
Species	Human
GeneID (Human)	4792
GeneID	4792
Immunogen	Synthetic phospho peptide corresponding to residues surrounding Y42 of human IkBa
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P25963
Protein Weight	36kDa
Shipping	Ice bag