


Phospho-Ikky (S376) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y6K9

Item Number: CM0001683



Introduction

Phospho-specific rabbit polyclonal antibody targeting IKKγ phosphorylated at Ser376, validated for Western blot and ELISA in human, mouse, and rat samples. Ideal for studying NF-κB signaling and immune response pathways.

[Learn More](#)

Product Name	Phospho-IKKγ-S376 Rabbit pAb
Remarks/Alias	IP; IP1; IP2; FIP3; IKKG; IPD2; NEMO; FIP-3; Fip3p; IMD33; SAIDX; AMCBX1; EDAID1; IKKAP1; ZC2HC9; IKK-gamma; Phospho-IKKγ-S376
Species	Human
Gene ID (Human)	8517
Gene ID	8517
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S376 of human IKKγ (NP_003630.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6K9
Protein Weight	48 kDa/56 kDa/37 kDa
Shipping	Ice bag