

Phospho-Insr-Y1189/Y1190 + Igf1R-Y1165/Y1166 Rabbit Pab For Wb, Elisa - P08069/P06213

Item Number: CM0000855



Introduction

Rabbit polyclonal antibody detecting dual phosphorylation sites Y1189/Y1190 on INSR and Y1165/Y1166 on IGF1R. Validated for WB and ELISA in human, mouse, rat samples. Ideal for studying insulin/IGF-1 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-INSR-Y1189/Y1190 + IGF1R-Y1165/Y1166 Rabbit pAb
Remarks/Alias	IGF1R/INSR; Phospho-INSR-Y1189/Y1190 + IGF1R-Y1165/Y1166
Species	Human
GeneID (Human)	3480/3643
GeneID	3480/3643
Immunogen	A synthetic phosphorylated peptide around Y1165 & Y1166 of human IGF1R/INSR (NP_000199.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08069/P06213
Protein Weight	154kDa/155kDa/156kDa
Shipping	Ice bag