


Phospho-Irak4-Ser346 Rabbit Mab For Wb, Elisa - Q9Nwz3

Item Number: CM0001395



Introduction

Phospho-specific rabbit monoclonal antibody against IRAK4 phosphorylated at Ser346. Validated for WB and ELISA, reacts with human IRAK4 (Q9NWZ3), a key kinase in innate immune signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRAK4-Ser346 Rabbit mAb
Remarks/Alias	IPD1; IMD67; REN64; IRAK-4; NY-REN-64
Species	Human
Gene ID (Human)	51135
Gene ID	51135
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Ser346 of human IRAK4 (NP_057207.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NWZ3
Protein Weight	52kDa
Shipping	Ice bag