


Phospho-Ire1-S724 Rabbit Pab - 075460

Item Number: CM0001877



Introduction

Rabbit polyclonal antibody targeting phospho-IRE1 (S724) for WB, IHC-P, ELISA applications. Validated in human, mouse, rat. Key tool for investigating ER stress and unfolded protein response pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRE1-S724 Rabbit pAb
Remarks/Alias	IRE1; IRE1P; IRE1a; hIRE1p; Phospho-IRE1-S724
Species	Human
GeneID (Human)	2081
GeneID	2081
Immunogen	A synthetic phosphorylated peptide around S724 of human IRE1alpha (NP_001424.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075460
Protein Weight	110 kDa
Shipping	Ice bag