

Phospho-Irf3-S386 Rabbit Pab For Wb And Elisa - Q14653

Item Number: CM0001275



Introduction

Rabbit polyclonal antibody against phosphorylated IRF3 at S386, validated for WB and ELISA applications. Reacts with human, mouse, and rat IRF3. Suitable for investigating interferon signaling and innate immune pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRF3-S386 Rabbit pAb
Remarks/Alias	IIAE7; Phospho-IRF3-S386
Species	Human
GeneID (human)	3661
GeneID	3661
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S386 of human IRF3 (NP_001562.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14653
Protein Weight	47kDa
Shipping	Ice bag