


Phospho-Irs1-S1101 Rabbit Pab For Wb, Elisa - P35568

Item Number: CM0000557



Introduction

Phospho-specific rabbit polyclonal antibody targeting IRS1 phosphorylated at Ser1101. Validated for Western blot and ELISA, cross-reacts with mouse. Ideal for insulin signaling and metabolic pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRS1-S1101 Rabbit pAb
Remarks/Alias	HIRS-1; Phospho-IRS1-S1101
Species	Human
Gene ID (Human)	3667
Gene ID	3667
Immunogen	A synthetic phosphorylated peptide around S1101 of human IRS1 (NP_005535.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P35568
Protein Weight	132kDa
Shipping	Ice bag